
Accolade® II

vástago de cadera

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Toolkit Return to menu

Scenario: Femoral Neck Resection

You are doing your first DSA case. The procedure has gone well so far. You have dissected and retracted the confluence of the piriformis and obturator internus tendons and performed the capsulotomy, but you run into difficulty during the osteotomy. While the hip may be dislocated, you still cannot adequately visualize if the resection has been made at the templated level.

Patient profile:
• 58-year-old female • BMI of 32

What should you do? Make the best selection.

A Decision A
Internally rotate the femur to 30 degrees.

B Decision B
Perform a subcapital osteotomy. Then replace the femoral straddlers to visualize the remaining neck and make a second neck cut. Assess the contained neck resections to the pre-op templated neck resection height.

C Decision C
Reduce the femur back into the acetabulum and make an in situ neck resection.

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
- National Center for Health Statistics. Health, United States, 2010: With Special Feature on Death and Dying, págs. 17 y 60.
- Kurtz S, Lau E, Ong K, Zhao K, Kelly MP, Bozic KJ. Future Young Patient Demand for Primary and Revision Joint Replacement. CORR 2009; 467:2606-2612.
- Haddad F. A Positive End to 2016. Bone Joint J. 2016;98-B:1569-70.
- Casper D, Kim G, Restrepo C, Parvizi J, Rothman R. Primary Total Hip Arthroplasty With an Uncemented Femoral Component. The Journal of Arthroplasty. Vol. 26, Edición 6, 2011; 838-841.
- McLaughlin JR, Lee KR. Total hip arthroplasty with an uncemented femoral component. J Bone Joint Surg 2008; 90:1290.
- Paul van der Voort, Bart G Pijls, Marc J Nieuwenhuijse, Jorrit Jasper, Marta Fiocco, Josepha W M Plevier, Saskia Middeldorp, Edward R Valstar & Rob G H H Nelissen (2015) Early subsidence of shape-closed hip arthroplasty stems is associated with late revision, Acta Orthopaedica, 86:5, 575-585, DOI: 10.3109/17453674.2015.1043832.
- Cooper J, et al. Distal Fixation of Proximally Coated Tapered Stems May Predispose to a Failure of Osteointegration. Journal of Arthroplasty. 2011. 10.1016/j.arth.2011.04.003.
- Abdel MP, et al. Epidemiology of periprosthetic fracture of the femur in 32 644 primary total hip arthroplasties. Bone Joint J 2016;98-B:461-7.
- Issa, K. et al. Radiographic Fit and Fill Analysis of a New Second-Generation Proximally Coated Cementless Stem Compared to its Predicate Design. Journal of Arthroplasty (2013). <http://dx.doi.org/10.1016/j.arth.2013.04.029>.
- Collopy, D. et al. A Prospective Roentgen Stereophotogrammetric Analysis (RSA) Study of the Stryker Accolade II Cementless Hipo Stem to 2 Years Post-Surgery. Resúmenes del AOA anual de 2016 Meeting-Cairns, QLD, 9-13 de octubre, 2016.
- Bono, J. et al. Total Hip Intraoperative Femur Fracture: Does a Tapered Wedge Design with a Size Specific Medial Curvature Reduce the Incidence of Intraoperative Femur Fracture? Presentación en la XLVII reunión anual de la EOA - Nueva Orleans, LA, 19-22 de octubre, 2016.
- Pierce T. et al. Second-Generation Versus First-Generation Cementless Tapered Wedge Femoral Stems. Orthopedics. 2015 38(9): 550-554.
- Kolisek, F. et al. Outcomes of Second-Generation Tapered Wedge Design. SurgTechnolInt. 4 de abril de 2016; XXVIII. pii: sti28/702.
- Presentación interna de Stryker 'Global Bone Morphology Study -Proximal Femur'. 10 de noviembre de 2016.
- Narzikul, BS. et al. Alteration in Geometry of Femoral Stem Results in Better Fit and Fill: Comparison of Accolade I vs. Accolade II. Presentación en la XLVII reunión anual de la EOA - Nueva Orleans, LA, 19-22 de octubre, 2016.
- Wuestemann T, Bastian A, Parvizi J, Nessler J, Kolisek F. A novel tapered hip stem design optimized for femoral fit in a wide array of bone types. EFFORT 2011.
- Boucher, F. et al. Preclinical Finite Element Analysis of a Novel Tapered Wedge Stem Optimized for Stability at a Reduced Length. Presentación de póster en EFFORT 2011.
- Lovell T, Hozack W, Kreuzer S, Merritt P, Nogler M, Puril, Wuestemann T, Bastian A. Influence of Stem Length on the Insertion Path in THR. Póster n.º 2040 en ORS 2009.
- Khanuja H, et al. Current Concepts Review: Cementless Femoral Fixation in Total Hip Arthroplasty. J Bone Joint Surg Am. 2011;93:500-9.
- Faizan A, et al. Development and Verification of a Cementless Novel Tapered Wedge Stem for Total Hip Arthroplasty, J Arthroplasty (2014), <http://dx.doi.org/10.1016/j.arth.2014.09.023>.
- Restrepo C, Parvizi J, Eslam Pour A, et al. Prospective Randomized Study of Two Surgical Approaches for Total Hip Arthroplasty. The Journal of Arthroplasty 2010; 25: 671-670.
- Vail T, Mariani E, Bourne M, Berger R, Meneghini M. Approaches in Primary Total Hip Arthroplasty. J Bone Joint Surgery 2009; 91:10.
- M.W. Zawadsky y cols./ Early Outcome Comparison Between the Direct Anterior Approach and the Mini-Incision Posterior Approach for Primary Total Hip Arthroplasty: 150 Consecutive Cases, The Journal of Arthroplasty 29 (2014) 1256-126.
- Datos de ventas internos (2017).

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Copyright © 2018 Stryker